

Application Note





Pipetting suspensions for process control in chemical production

Anders Karlsson 1, Dag Dacke 2, Renzo Picenoni 3, Guenter Boehm 3



¹ AkzoNobel Pulp and Performance Chemicals AB, Sundsvall, Sweden

² Kovalent AB, Västra Frölunda Sweden

³ CTC Analytics AG, Zwingen, Switzerland

Introduction

AkzoNoble is the producer of Expancel. Expancel microspheres (fig. 1) are small spherical thermoplastic particles. The microspheres consist of a polymer shell encapsulating a gas. When heated, the internal pressure from the gas increases and the thermoplastic shell softens, resulting in a dramatic increase of the volume of the microspheres. The gas remains inside the spheres.



Fig 1: Enlargement of a heated Expancel microspehere

Process control

During the production of Expancel samples are regularly analysed for a number of critical parameters (fig. 3). So far this sample preparation has been a manual process. To increase productivity a semi-automated sample preparation has been evaluated. 100 mL samples are still taken manually at line. However, the processing and analysis of the samples should then be fully automated.

A critical point in the evaluation of the automated sample prep was to find out if pipetting can be used to reliably transfer defined amounts of the sample suspension (see cover photo).





Fig. 3: Process steps for the at-line analysis of Expancel production batches



Fig. 4: Out-line of hardware set up for automation



Fig.5 PAL RTC setup for automated weighing of pipetted samples



Application Note: Weighing samples automatically: Gravimetric dilutions of highly viscous samples

| Weigh | ht of transferred suspension [mg] | | _ | | |
|---------|---------------------------------------|---------|---------------------|----------------------|--------------------|
| 18.5000 | | Sample | Weight Tara [mg] | Weight Gross [mg] | Weight Net [mg] |
| | | | 1 17.64339 | 17.65998 | 16.5900 |
| 18.0000 | | | 2 17.54938 | 17.5668 | 17.4200 |
| | | | 3 17.67467 | 17.69065 | 15.9800 |
| 17.5000 | • • • • • • • • • • • • • • • • • • • | | 4 17.58607 | 17.60299 | 16.9200 |
| | | | 5 17.60297 | 17.62001 | 17.0400 |
| 17.0000 | | | 6 17.61459 | 17.63162 | 17.0300 |
| | • • | | 7 17.69568 | 17.71218 | 16.5000 |
| 16.5000 | | | 8 17.53305 | 17.54977 | 16.7200 |
| | • | | 9 17.45578 | 17.4725 | 16.7200 |
| 16.0000 | | | 10 17.53061 | 17.54695 | 16.3400 |
| | | | 11 17.70394 | 17.72066 | 16.7200 |
| 15 5000 | | | 12 17.70634 | 17.72363 | 17.2900 |
| 10.0000 | | Average | | | 16.77250 |
| 15.0000 | | SD | | | 0.40200 |
| | | %RSD | | | 2.40 |
| | | | | | |
| 14.5000 | | | | | |
| | | | | | |
| 14.0000 | | | | | |
| 1 | 1 3 5 7 9 11 | 13 | | | |

Fig 6: Typical results for a particular production batch

PAL Pipette Tool – Essential for sample preparation

Pipetting liquids with tips is a common practice in many applications. The most important feature is probably the fact that carryover is zero by design because the used tip is discarded after every sample transfer. Moreover, functionalized tips (e.g. for desalting or protein capture) allow for effective sample work-up.

The new PAL Pipette Tool significantly enhances the PAL System sample prep platform. In combination with the complete range of tools and modules the PAL Pipette Tool facilitates seamless lab workflows, especially in combination with the PAL DeCapper Module:



Only the PAL RTC can switch automatically between pipettes and positive displacement syringes. Pipettes are free from carryover by design. Chemically inert glass syringes are best when highest accuracy and compatibility with organic solvents is required.



Application areas:

- Serial dilutions
- Micro SPE in a tip
- Protein digestion in a tip
- Reformatting of samples
- Pipetting of suspensions
 - (e.g. magnetic beads) and emulsions

Technical data, ordering information:

0-200 μL P/N PAL3-TH-Pip-200 0-1000 μL P/N PAL3-TH-Pip-1000

PAL3 Firmware: 2.4.x or later

Detailed specifications are available.



Fig. 7: Wash module and adapter for LC-injection port with Pipette Tool (PN PAL3-WashPip)

Highlights:

- Carryover free sample handling
- Complete workflows in combination with PAL Tool and Modules
- Direct injection from the Pipette Tool into HPLC valves (see fig. 1)
- Glass syringe (positive displacement) and Pipette Tool on one system
- Tested PAL Consumables like pipette tips and racks enhance the performance further



Fig. 8: Pipette Tool placed in Park Station



Pipette Tool, Gravimetric Specifications

Repeatability, accuracy and linearity, single dispense, 200 μL Pipette Tool



Repeatability, accuracy and linearity, single dispense, 1000 μL Pipette Tool

n=7



| Dispense Volume | Repeatability | Accuracy | Linearity |
|-----------------|------------------|-------------|-----------------|
| 20 µL | $RSD \le 1.5 \%$ | ≤2 % | |
| 100 µL | $RSD \le 0.5~\%$ | ≤1% | $R^2 \ge 0.999$ |
| 200 µL | $RSD \le 0.5~\%$ | $\leq 1 \%$ | |

| Dispense Volume | Repeatability | Accuracy | Linearity |
|-----------------|------------------|-----------------|-----------------|
| 100 µL | $RSD \le 0.5~\%$ | ≤ 0.5 % | |
| 500 µL | $RSD \le 0.3~\%$ | ≤ 0.5 % | $R^2 \ge 0.999$ |
| 1000 μL | $RSD \le 0.3 \%$ | $\leq 0.5 \ \%$ | |

Conditions:

Pipette Tool 200 μ L, measured with two different tools Liquid: water, 23 °C, density 0.997538 g/mL Balance used: Precisa, XB 120A (d=0.1 mg) PAL System 200 μ L tips

| Liquid Class Parameter | Liquid Class: Water (Single Dispense) |
|---------------------------|---------------------------------------|
| Airgap, front [µL] | 0 |
| Airgap, rear [µL] | 100.0 |
| Aspirate delay [ms] | 300 |
| Aspirate flow rate [µL/s] | 100.0 |
| Correction factor | 1.04 |
| Dispense delay [ms] | 100 |
| Dispense flow rate [µL/s] | 100.0 |
| Volume offset [µL] | 1.6 |
| Volume, front [µL] | 0 |
| Volume, rear [µL] | 0 |

PAL Sample Control Method: Aspirate Dispense with Liquid Class_Rev1

Conditions:

Pipette Tool 1000 μ L, measured with two different tools Liquid: water, 23 °C, density 0.997538 g/mL Balance used: Precisa, XB 120A (d=0.1 mg) PAL System 1000 μ L tips

| Liquid Class Parameter | Liquid Class: Water (Single Dispense) |
|---------------------------|---------------------------------------|
| Airgap, front [µL] | 0 |
| Airgap, rear [µL] | 20.0 |
| Aspirate delay [ms] | 300 |
| Aspirate flow rate [µL/s] | 100.0 |
| Correction factor | 1.026 |
| Dispense delay [ms] | 100 |
| Dispense flow rate [µL/s] | 500.0 |
| Volume offset [µL] | 5.9 |
| Volume, front [µL] | 0 |
| Volume, rear [µL] | 0 |

PAL Sample Control Method: Aspirate Dispense with Liquid Class_Rev1

Conclusions:

- The PAL RTC with Pipette Tool allowed for the precise pipetting of the Expancel suspensions.
- The precision of the pipetting step has been assessed by automatic weighing.
- Batch to batch variations showed RSDs 1 4%.
- The automation of the sample preparation allows for a significant increase in productivity since sample preparation during the night or over weekends is now possible when no skilled lab personnel is available.

Legal Statements

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described in this document at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall CTC Analytics AG be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.

© 2018 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Imprint

Date of print: 04.2018

CTC Analytics AG Industriestrasse 20 CH-4222 Zwingen Switzerland T +41 61 765 81 00 Contact: info@ctc.ch