

LOW LEVEL ANALYSIS OF PHENOLS BY GC/MSMS AFTER ON-LINE ACETYLATION WITH A ROBOTIC SAMPLE HANDLER

Cedric WISSEL, Aldo LUNARDI and Joeri VERCAMMEN

Phenols comprise a class of chemical components consisting of a hydroxyl moiety, which is directly bound to an aromatic ring, and a range of alkyl, chlorine and nitro-substituents. Phenols are applied in various chemical processes and their presence in the environment needs to be monitored carefully. Analysis is generally carried out by means of GC/MS after acetylation to prevent severe peak tailing that compromises chromatographic separation from occurring. Unfortunately, sample preparation is tedious and time-consuming and highly prone to user-induced variability.

In this report we present the performance of a fully automated method in which phenols are acetylated according to standard reference procedures and analyzed using GC/MSMS after in-vial liquid/liquid extraction (LLE) and large volume PTV injection.

MATERIALS AND METHODS

Target components:

Table 1. Overview of target components.

| Component | Component |
|-------------------------|---------------------------|
| 2-methylphenol | 4-chloro-3-methylphenol |
| 3-methylphenol | 2,5-dichlorophenol |
| 4-methylphenol | 2,4-dichlorophenol |
| 2-chlorophenol | 3,5-dichlorophenol |
| 2,6-dimethylphenol | 2,3,5-trimethylphenol |
| 2-ethylphenol | 2,3-dichlorophenol |
| 3-chlorophenol | 3,4-dichlorophenol |
| 2,5-dimethylphenol | 2,4,6-trichlorophenol |
| 4-chlorophenol | 2,3,6-trichlorophenol |
| 2,4-dimethylphenol | 2,3,5-trichlorophenol |
| 3-ethylphenol | 2,4,5-trichlorophenol |
| 2-isopropylphenol | 2,3,4-trichlorophenol |
| 2,3-dimethylphenol | 3,4,5-trichlorophenol |
| 4-ethylphenol | 2,3,5,6-tetrachlorophenol |
| 3,5-dimethylphenol | 2,3,4,6-tetrachlorophenol |
| 3,4-dimethylphenol | 2,3,4,5-tetrachlorophenol |
| 4-chloro-2-methylphenol | pentachlorophenol |
| 2,6-dichlorophenol | bisfenol-A |

Standards:

Stock solutions of native and labeled components were provided by the customer. Standards were prepared by adding the native components to drinking water before extraction/derivatization. An overview of standard concentrations is given in Table 2. Please note that at this stage only one labeled component was applied, i.e. 3,4-dichlorophenol-D3.

Table 2. Overview of standards and concentrations.

| Reference | Concentration (µg/L) |
|---------------------------|----------------------|
| Standard 1 | 0.01 |
| Standard 2 | 0.02 |
| Standard 3 | 0.05 |
| Standard 4 | 0.10 |
| Standard 5 | 0.20 |
| Standard 6 | 0.50 |
| Standard 7 | 1.00 |
| Drinking water addition 1 | 0.10 |
| Drinking water addition 2 | 0.10 |
| Drinking water addition 3 | 0.10 |
| Drinking water addition 4 | 0.10 |
| Drinking water addition 5 | 0.10 |
| Surface water addition | 0.10 |
| Ground water addition | 0.10 |

Autosampler settings:

In-vial extraction and derivatization was carried out by means of a Thermo TriPlus RSH robotic sampler, which is placed on top of the GC/MSMS instrument (Thermo Trace 1300 GC and TSQ 8000 Evo, see Figure 1).

The device is equipped with the following (modular) tools:

- Tool stations for 10 mL, 1 mL and 100 µL tools
- Vortex option
- Fast Wash Station
- Solvent reservoirs for 100 mL vials
- Solvent reservoir for 10 mL vials
- Trays for 20 mL and 2 mL vials

The software control is completely embedded in the Xcalibur/TraceFinder instrument method. Method overlapping is enabled; the next sample is prepared during the analysis of the previous. Instrumental settings are described in Table 3.



Figure 1. Instrumental set-up.

RESULTS AND DISCUSSION




Table 3. Validation results.

| Compound | %RSD (n=5) | R ² | Surface water, spike @ 0.1 µg/L |
|---------------------------|------------|----------------|---------------------------------|
| 2-methylphenol | 11.4 | 0.9999 | 0.11 |
| 3-methylphenol | 12.9 | 0.9987 | 0.10 |
| 4-methylphenol | 9.5 | 0.9998 | 0.11 |
| 2-chlorophenol | 7.4 | 0.9998 | 0.11 |
| 2,6-dimethylphenol | 7.1 | 0.9990 | 0.12 |
| 2-ethylphenol | 6.7 | 0.9998 | 0.11 |
| 3-chlorophenol | 8.4 | 0.9999 | 0.11 |
| 2,5-dimethylphenol | 7.3 | 0.9994 | 0.11 |
| 4-chlorophenol | 9.3 | 0.9998 | 0.11 |
| 2,4-dimethylphenol | 6.3 | 0.9997 | 0.12 |
| 3-ethylphenol | 7.1 | 0.9999 | 0.11 |
| 2,3-dimethylphenol | 7.0 | 0.9997 | 0.11 |
| 3,5-dimethylphenol | 8.1 | 0.9986 | 0.11 |
| 4-ethylphenol | 8.1 | 0.9986 | 0.11 |
| 3,4-dimethylphenol | 7.6 | 0.9997 | 0.11 |
| 4-chloro-2-methylphenol | 7.2 | 0.9997 | 0.11 |
| 2,6-dichlorophenol | 6.5 | 0.9999 | 0.11 |
| 4-chloro-3-methylphenol | 6.9 | 0.9998 | 0.11 |
| 2,5-dichlorophenol | 7.0 | 0.9999 | 0.11 |
| 2,4-dichlorophenol | 6.5 | 1.0000 | 0.10 |
| 3,5-dichlorophenol | 6.7 | 1.0000 | 0.11 |
| 2,3-dichlorophenol | 7.8 | 0.9999 | 0.10 |
| 3,4-dichlorophenol | 7.2 | 0.9997 | 0.10 |
| 2,4,6-trichlorophenol | 6.6 | 1.0000 | 0.11 |
| 2,3,6-trichlorophenol | 6.1 | 0.9999 | 0.11 |
| 2,3,5-trichlorophenol | 6.4 | 0.9999 | 0.10 |
| 2,4,5-trichlorophenol | 6.4 | 1.0000 | 0.09 |
| 2,3,4-trichlorophenol | 7.2 | 0.9999 | 0.10 |
| 3,4,5-trichlorophenol | 7.2 | 1.0000 | 0.10 |
| 2,3,5,6-tetrachlorophenol | 7.1 | 0.9999 | 0.11 |
| 2,3,4,6-tetrachlorophenol | 7.3 | 1.0000 | 0.10 |
| 2,3,4,5-tetrachlorophenol | 8.5 | 0.9999 | 0.10 |
| Pentachlorophenol | 9.1 | 0.9997 | 0.11 |

CONCLUSIONS

A fully automated large volume GC/MSMS method was developed for the analysis of residual phenols in water after on-line derivatization with acetic anhydride. The method was fully validated with linearity coefficients > 0.998 for all components by using only one labelled internal standard. Method precision was < 10% for all components with minor deviation for the most volatile phenols in the target list. An additional, volatile, ISTD will act as a more effective control and reduce this phenomenon. Recovery was tested for two matrices with results fulfilling typical laboratory requirements.