



Automated In-Vial Extraction of Surface Waters for the Analysis of Mineral Oil by GC-FID

Mathias Vanlancker
Analytical support and method development

Sample Prep

Manual Procedure


- Add 12 mL of surface water to a 20 mL headspace vial
- Place on the Triplus RSH robotic sampler
- Prep 2ml vials with a scoop of Florisil

Automated Procedure

- 3 (or 4) mL of hexane is added
- C10 and C40 ISTD markers are added by the RSH
- Vortex
- Inject 50 μ L onto GC-FID (LV-PTV injection)
- Transfer 1.5 mL to a 2 mL vial containing a scoop of Florisil (can be done for positives only or all samples)
- Vortex
- Inject 50 μ L onto the GC-FID

Calibration Curve

A calibration curve is made between 0.0625 ppm and 2.5 ppm (3 mL extraction)
Or a calibration curve between 0.0833 ppm and 3.33 ppm (4 mL extraction)



$R^2 = 0.9994$

No offset!

Calibration Curve

	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6
Rep 1	0.046	0.116	0.239	0.620	1.274	2.509
Rep 2	0.045	0.118	0.230	0.596	1.239	2.473
Rep 3	0.045	0.113	0.230	0.612	1.222	2.529
Average (ppm)	0.046	0.115	0.233	0.609	1.245	2.504
St. Dev. (ppm)	0.00069	0.002537	0.005104	0.012302	0.026522	0.028784
RSD (%)	1.5	2.2	2.2	2.0	2.1	1.1
LOD (ppm)	0.002					
LOQ (ppm)	0.007					

LOD = 3 * St. Dev.


LOQ = 10 * St. Dev.

Residuals

	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6
Rep 1	-25.72	-7.4	-4.37	-0.77	1.93	0.37
Rep 2	-27.8	-6	-7.85	-4.63	-0.84	-1.09
Rep 3	-27.4	-10	-7.96	-2.03	-2.22	1.17

Overlay chromatogram

Benelux RIVM Hydrocarbon Standard



Method Summary

- Cycle time = 21 min! -> 68 analysis/24h
- No pre-column, 15 m column direct from PTV to FID
- Only 3 or 4 mL Hexane per sample
- RSD < 3%
- $R^2 > 0.999$
- Low LOD 0.002ppm
- Low LOQ 0.007ppm



Conclusions

- Fast method with Prep Ahead function
- PTV analysis: Simple maintenance and more robust against matrix compared to Oncolumn
- Low solvent use
- Excellent initial validations
- No SPE necessary
- Only 1 software: **Chromeleon**
- Completely automated due to smart sequencing



Add some Q
to your lab

inge.de.dobbeleer@interscience.nl
m.vanlancker@interscience.be

