


SampleQ: Laboratory automation and more



New: Exetainer tray



Exetainer screw capped vials from Labco Ltd. are used routinely for breath analysis in gastroenterology clinics, pharmaceutical laboratories and research test centers. They are highly valued for ^{13}C Helicobacter pylori Testing, Breath Test Kits, Gas and Liquid Testing, Specimen collection, storage and much more. Although manufactured to fit most analyzers, they were incompatible with GC and ICP-MS autosamplers until recently.

The Exetainer Tray Holder is a new option available from SampleQ that can be installed on the Thermo Scientific TriPlus RSH autosampler and the CTC PAL-3 models to sample directly from 12 mL Exetainer vials. It consists of a fixed tray that holds 96 samples and a vial retaining plate that can be installed without any tool.

Features include:

- High capacity tray
- Detailed installation & user manual
- Secure locking of vials with an easy to use retaining plate

More information on
www.sampleq.com

