

# SampleQ: Laboratory automation and more



## New: ZSC module



SampleQ has developed a new option for the RSH/RTC PAL-Autosampler that allows **zone-specific cooling** of a single tray cassette. Contrary to the classic solution, in which vials are stored in a bulky stack configuration, our approach is based on a standard triple cassette trayholder.

Benefits include:

- Compact design
- Eco-friendly
- Condensation-free operation due to active ventilation
- Easy access to all sample positions
- Variable temperature zones per tray holder

More information on  
[www.sampleq.com](http://www.sampleq.com)

