


SampleQ: Laboratory automation and more



New: Large Vial Option



SampleQ has developed a new option for the RTC PAL Auto-sampler that allows storage and transport of larger vials.

Now the 40 mL or VOA-type of vials you are using in your daily lab routine are easily transferred to this powerful automated sample prep platforms.

Features and benefits:

- Up to 36 samples storage rack
- Suitable for 40 mL or VOA vials
- Automated vial detection checks samples are present prior to analysis
- Standard: magnetic transport (requires metal caps)
- Optional: Mechanical transport with gripper tool
- No sample transfer needed to switch to the PAL RTC platform
- Reduces human error
- Cost efficient use of vials

More information on
www.sampleq.com

